

» Verification Cabinet for Lead Block



- » The light source can be accommodated freely, and the center vertical point can be found conveniently
- » The distance between the light source and pallet can be adjusted up and down freely
- » The pallet can be regulated according to the size of the blocking trays
- » The guiding rail can slide the blocks more freely even carrying 100kg lead

» Manual Hot-wire Cutter




- » Used for x-ray treatment
- » The maximum distance between light source and film is 160cm
- » Balance block on both sides are used for balance adjustment
- » here are three different lengths of drawing arms
- » The adjustable device on the bottom can be adapted the uneven groundv

» Shoulder Retractor



» Water Bath



- » Stainless-steel material:
- » Shape size: 650*600*200mm