


P M

Patient Marker



1,1.2,1.8,2,2.5 MM,DIA



Blue CT Markers



Cross Markers



Multimodality markers



Line Marker



CT Mark Wire, 2.0mm Wire



Sure Mark



Vaginal Dilator



Gold fiducial markers

STRAIGHT & CURVED BODY CALIPER



The Rectal Marker is a flexible tube packed with lead balls spaced at 1 cm intervals from center to center. The CT Rectal Marker has Teflon balls with same spacing. The Rectal Marker is used to accurately obtain both the rectum position and magnification by counting the balls. This determines the rectum location relative to the radiation field. An adjustable anus marker can be utilized during simulation or port filming. The Rectal Marker can be used in the vagina or as an external marker