

Block Casting



1. Die steel material;
2. It has screws for horizontal adjustment.

Radiotherapy Bower



YC-HM176

1. Used for x-ray treatment;
2. Reflection and scattering screen at both front and back;
3. Two hand grips at both sides;
4. Casters at the bottom are easy to mobile.

Horizontal Radiotherapy Bed



1. Used for electronic-ray treatment;
2. Steel material bed and aluminium alloy material light box;
3. Reflection and scattering screen at both front and back, and the thickness of the screens is 10mm;
4. The height range of the treatment bed is 176-216mm.