

# Laser

## » Fixed & Movabel Laser



MICRO™ DIODE  
LASER



### Red Fixed Laser Systems

Diacor offers three configurations for our Centralite red diode lasers that address common treatment requirements and make ordering from Diacor fast and easy.



Red Fixed Laser LS-1  
Four fixed laser system



Red Fixed Laser LS-2  
Three fixed laser system.




Red Fixed Laser LS-3  
Three red fixed laser system with overhead mount.



### Green Fixed Laser System


Green lasers deliver the same accuracy as red diode lasers while projecting an image that is easier to see across a variety of skin tones.



Green Fixed Laser LS-4  
Four green laser system that offers vertical and horizontal mounting options.



Green Fixed Laser LS-5  
Three green fixed laser system offers horizontal and vertical mounting options.



Green Fixed Laser LS-6  
Three fixed green laser system that offers two cross lasers and ceiling mount.

## » Technical Parameter

Each laser's parameters are as following:

- » Laser diode
- » Line Width:  $\leq 1\text{mm}$  (in the distance 3~5m)
- » Line Length  $\geq 2\text{m}$  (at 3 m distance)
- » Wavelength: 635nm