

S

P

P

» Stamps for Radiotherapy



- » 40 pieces in per set
- » Including all body organ shapes
- » Packed in Aluminum Box

» Electrical Convection Oven

For softening thermoplastic material

- » Heat circulation, convection principal
- » Touch Display for temperature adjustment
- » Simple...Easy...No Preheating
- » No more wet towels and steamy treatment rooms
- » Cleaner and safer than waterbaths
- » Reduced risk of bacterial contamination

- » Space-saving
- » Compact design
- » Non-stick surface
- » Insulated case
- » Dry heating eliminates the inconvenience of using water
- » Digital temperature setting in Celsius or Fahrenheit
- » Auto shut-down after 2 hours—Turns off after 6 hours.
- » Finish, White baked enamel (on cold rolled steel),stainless steel trim Outer dimensions: 31x23x5.5 in Weight:46 lb (21 Kg) Voltage:110V and 220V

